


BIOL647
Digital Biology

Rodolfo Aramayo



FASTQ Format

Ben Langmead




Department of Computer Science

Aaron Quinlan
Departments of Human Genetics and Biomedical Informatics
USTAR Center for Genetic Discovery
University of Utah

FASTQ Format

The majority of DNA sequencing technologies produce a FASTQ file



Name

Sequence (ignore)

Base qualities

@ERR194146.1 HSQ1008:141:D0CC8ACXX:3:1308:20201:36071/1

ACATCTGGTTCCTACTTCAGGGCCATAAAGCCTAAATAGCCCACACGTTCCCCTTAAAT

+

?@@FFBFFDDHHBCEAFGEGIIDHGH@GDHHHGEHID@C?GGDG@FHIGGH@FHBEG:G

FASTQ Format

A "standard" format for storing and defining sequences from next-generation sequencing technologies http://en.wikipedia.org/wiki/FASTQ_format



Name

Sequence (ignore)

Base qualities

@ERR194146.1 HSQ1008:141:D0CC8ACXX:3:1308:20201:36071/1

ACATCTGGTTCCTACTTCAGGGCCATAAAGCCTAAATAGCCCACACGTTCCCCTTAAAT

+

?@@FFBFFDDHHBCEAFGEGIIDHGH@GDHHHGEHID@C?GGDG@FHIGGH@FHBEG:G

FASTQ Format

The FASTQ format's sequence identifier (first line of each record)

Old format

@HWUSI-EAS100R:6:73:941:1973#0/1

| HWUSI-EAS100R | the unique instrument name | | | | | |
|---------------|---|--|--|--|--|--|
| 6 | flowcell lane | | | | | |
| 73 | tile number within the flowcell lane | | | | | |
| 941 | 'x'-coordinate of the cluster within the tile | | | | | |
| 1973 | 'y'-coordinate of the cluster within the tile | | | | | |
| #0 | #0 index number for a multiplexed sample (0 for no indexing) | | | | | |
| /1 | the member of a pair, /1 or /2 (paired-end or mate-pair reads only) | | | | | |

New format

@EAS139:136:FC706VJ:2:2104:15343:197393 1:Y:18:ATCACG

| EAS139 | the unique instrument name | | | |
|---------|--|--|--|--|
| 136 | the run id | | | |
| FC706VJ | the flowcell id | | | |
| 2 | flowcell lane | | | |
| 2104 | tile number within the flowcell lane | | | |
| 15343 | 'x'-coordinate of the cluster within the tile | | | |
| 197393 | 'y'-coordinate of the cluster within the tile | | | |
| 1 | the member of a pair, 1 or 2 (paired-end or mate-pair reads only) | | | |
| Y | Y if the read is filtered, N otherwise | | | |
| 18 | 0 when none of the control bits are on, otherwise it is an even number | | | |
| ATCACG | index sequence | | | |
| | | | | |

FASTQ Format

FASTQ quality scores: estimate of confidence in each base (sequencing technologies make errors!)

```
@SEQ_ID
GATTTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGTTT
+
!''*((((***+))%%++)(%%%).1***-+*''))**55CCF>>>>>CCCCCCC65
```

Qualities are based on the Phred scale and are encoded:



$$Q = -10*log10(Perr)$$

Note:

The Ph in Phred comes from Phil Green, the inventor of the encoding

FASTQ Format


Phred quality score calculation

$$Q = -10*log10(Perr)$$

| Error probability | log10(Per) | Phred quality score | | | |
|-------------------|------------|---------------------|--|--|--|
| 1 | 0 | 0 | | | |
| 0.1 | -1 | 10 | | | |
| 0.01 | -2 | 20 | | | |
| 0.001 | -3 | 30 | | | |
| 0.0001 | -4 | 40 | | | |

FASTQ Format

A higher quality score is better (>=20 is considered "good")



FASTQ Format

Historically, FASTQ has had different encoding schemes for encoding PHRED quality scores

```
.......
   !"#$%&'()*+,-./0123456789:;<=>?@ABCDEFGHIJKLMNOPQRSTUVWXYZ[\]^_`abcdefghijklmnopqrstuvwxyz{|}~
33
                                  104
                                             126
0.2.....41
                                          Current encoding:
S - Sanger Phred+33, raw reads typically (0, 40)
        Solexa+64, raw reads typically (-5, 40)
X - Solexa
                                          ! = quality 0
I - Illumina 1.3+ Phred+64, raw reads typically (0, 40)
J - Illumina 1.5+ Phred+64, raw reads typically (3, 40)
                                          J = quality 41
 with 0=unused, 1=unused, 2=Read Segment Quality Control Indicator (bold)
  (Note: See discussion above).
L - Illumina 1.8+ Phred+33, raw reads typically (0, 41)
```

FASTQ Format


Quality score encoding based on ASCII table


| ١ | Dec | Hex | Char | Dec | Hex | Char | Dec | Hex | Char | Dec | Hex | Char |
|---|-----|-----|------------------|-----|-----|-------|-----|-----|------|-----|-----|------|
| H | 0 | 00 | Null | 32 | 20 | Space | 64 | 40 | 0 | 96 | 60 | ` |
| | 1 | 01 | Start of heading | 33 | 21 | ļ. | 65 | 41 | A | 97 | 61 | a |
| | 2 | 02 | Start of text | 34 | 22 | " | 66 | 42 | в | 98 | 62 | b |
| | 3 | 03 | End of text | 35 | 23 | # | 67 | 43 | С | 99 | 63 | c |
| | 4 | 04 | End of transmit | 36 | 24 | Ş | 68 | 44 | D | 100 | 64 | d |
| | 5 | 05 | Enquiry | 37 | 25 | 26 | 69 | 45 | E | 101 | 65 | e |
| | 6 | 06 | Acknowledge | 38 | 26 | ٤ | 70 | 46 | F | 102 | 66 | £ |
| | 7 | 07 | Audible bell | 39 | 27 | 1 | 71 | 47 | G | 103 | 67 | g |
| | 8 | 08 | Backspace | 40 | 28 | (| 72 | 48 | H | 104 | 68 | h |
| | 9 | 09 | Horizontal tab | 41 | 29 |) | 73 | 49 | I | 105 | 69 | i |
| | 10 | OA | Line feed | 42 | 2A | * | 74 | 4A | J | 106 | 6A | ژ |
| | 11 | ов | Vertical tab | 43 | 2B | + | 75 | 4B | K | 107 | 6B | k |
| | 12 | OC. | Form feed | 44 | 2 C | , | 76 | 4C | L | 108 | 6C | 1 |
| | 13 | OD | Carriage return | 45 | 2 D | - | 77 | 4D | M | 109 | 6D | m |
| | 14 | OE | Shift out | 46 | 2 E | | 78 | 4E | N | 110 | 6E | n |
| | 15 | OF | Shift in | 47 | 2 F | / | 79 | 4F | 0 | 111 | 6F | 0 |
| | 16 | 10 | Data link escape | 48 | 30 | 0 | 80 | 50 | P | 112 | 70 | p |
| | 17 | 11 | Device control 1 | 49 | 31 | 1 | 81 | 51 | Q | 113 | 71 | q |
| | 18 | 12 | Device control 2 | 50 | 32 | 2 | 82 | 52 | R | 114 | 72 | r |
| | 19 | 13 | Device control 3 | 51 | 33 | 3 | 83 | 53 | ສ | 115 | 73 | s |
| | 20 | 14 | Device control 4 | 52 | 34 | 4 | 84 | 54 | Т | 116 | 74 | t |
| | 21 | 15 | Neg. acknowledge | 53 | 35 | 5 | 85 | 55 | U | 117 | 75 | u |
| | 22 | 16 | Synchronous idle | 54 | 36 | 6 | 86 | 56 | v | 118 | 76 | v |
| | 23 | 17 | End trans, block | 55 | 37 | 7 | 87 | 57 | ឃ | 119 | 77 | w |
| | 24 | 18 | Cancel | 56 | 38 | 8 | 88 | 58 | x | 120 | 78 | x |
| | 25 | 19 | End of medium | 57 | 39 | 9 | 89 | 59 | Y | 121 | 79 | У |
| | 26 | 1A | Substitution | 58 | ЗА | : | 90 | 5A | Z | 122 | 7A | z |
| | 27 | 1B | Escape | 59 | зв | ; | 91 | 5B | [| 123 | 7В | { |
| | 28 | 1C | File separator | 60 | 3 C | < | 92 | 5C | ١ | 124 | 7C | I |
| | 29 | 1D | Group separator | 61 | ЗD | = | 93 | 5D |] | 125 | 7D | } |
| | 30 | 1E | Record separator | 62 | ЗE | > | 94 | 5E | ٨ | 126 | 7E | ~ |
| | 31 | 1F | Unit separator | 63 | ЗF | ? | 95 | 5F | | 127 | 7F | |

Formula for getting PHRED quality from encoded quality:

$$Q = (ascii(char) - 33)$$

Example:





BIOL647
Digital Biology

Rodolfo Aramayo

